Date prepared: May 28, 2012

# **Material Safety Data Sheet**

Date revised: N/A (first edition)

Section 1 - General Information		
Chemical Name:	Chlortetracycline hydrochloride	
CAS Reg. Number:	64-72-2	
Catalog Number:	1ST4027	
Manufacturer's Name:	A Chemtek, Inc.	
	100 Barber Avenue	
	Worcester, MA 01606, USA	
<b>Telephone Number:</b>	508-471-4121	
Fax Number:	508-845-9201	

For R & D use only

Section 2 - Hazardous Ingredients/Identity Information

**Emergency Overview OSHA** Hazards Target Organ Effect, Irritant, Teratogen **Target Organs** Liver, Teeth, Bone **GHS Classification** Acute toxicity, Oral (Category 5) Skin irritation (Category 2) Eye irritation (Category 2A) Specific target organ toxicity - single exposure (Category 3) GHS Label elements, including precautionary statements

Pictogram	$\mathbf{V}$	
Signal word	Warning	
Hazard statement(s)		
H303	May be harmful if swallowed.	
H315	Causes skin irritation.	
H319	Causes serious eye irritation.	
H335	May cause respiratory irritation.	
Precautionary stateme	ent(s)	
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.	
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if	
	present and easy to do. Continue rinsing.	
HMIS Classification		
Health hazard:	2	
Chronic Health H	Hazard: *	
Flammability:	0	
Physical hazards	: 0	
NFPA Rating		
Health hazard:	2	
Fire:	0	
Reactivity Hazar	d: 0	
Potential Health Eff	ects	
Inhalation	May be harmful if inhaled. Causes respiratory tract irritation.	
Skin	May be harmful if absorbed through skin. May cause skin irritation.	
Eyes	May cause eye irritation.	
Ingestion	May be harmful if swallowed.	

# Section 3 - Composition/Information on Ingredients

Formula:	C <sub>22</sub> H <sub>23</sub> ClN <sub>2</sub> O <sub>8</sub> ·HCl
Molecular Weight:	515.34g/mol
CAS-No.	64-72-2

**Section 4 - First Aid Measures** If inhaled

# MSDS

Date prepared: May 28, 2012

# **Material Safety Data Sheet**

Date revised: N/A (first edition)

If breathed in, move person into fresh air. If not breathing give artificial respiration. In case of skin contact Wash off with soap and plenty of water. In case of eye contact Flush eyes with water as a precaution. If swallowed Never give anything by mouth to an unconscious person. Rinse mouth with water.

# **Section 5 - Fire-Fighting Measures**

Suitable extinguishing media
Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.
Special protective equipment for fire-fighters
Wear self contained breathing apparatus for fire fighting if necessary.
Hazardous combustion products
Hazardous decomposition products formed under fire conditions. - Carbon oxides, nitrogen oxides (NOx), Hydrogen chloride gas, Hydrogen fluoride

# Section 6 - Accidental Release Measures

Personal precautions	Avoid dust formation.
Environmental precautions	Do not let product enter drains.
Methods for cleaning up	Sweep up and shovel. Keep in suitable, closed containers for disposal.

# Section 7 - Handling and Storage

HandlingProvide appropriate exhaust ventilation at places where dust is formed. Normal measures for<br/>preventive fire protection.StorageKeep container tightly closed in a dry and well-ventilated place.<br/>Recommended storage temperature: 2 - 8 °C

# Section 8 - Exposure Controls/Personal Protection

Contains no substances with occupational exposure limit values.

### Personal protective equipment

#### Eye/face protection

Safety glasses with side-shields conforming to EN166. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

# Skin protection

Complete suit protecting against chemicals. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

#### **Body Protection**

The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

#### **Respiratory protection**

For nuisance exposures use type P95 (US) or type P1 (EU EN 143) particle respirator.For higher level protection use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

# Section 9 - Physical and Chemical Properties

Solid
not available
Melting point/range: 243 - 247 °C (469 - 477 °F)
not available

**MSDS** 

Date prepared: May 28, 2012

# **Material Safety Data Sheet**

Date revised: N/A (first edition)

Other safety information not available

 Section 10 - Stability and Reactivity

 Storage stability
 Stable under recommended storage conditions

 Materials to avoid
 Strong oxidizing agents

 Conditions to avoid
 not available

 Hazardous decomposition products
 Image: Strong oxidizing agents

 Hazardous decomposition products
 Formed under fire conditions. - Carbon oxides, nitrogen oxides (NOx),

 Hydrogen chloride gas
 Other decomposition products - no data available

# **Section 11 - Toxicological Information**

Acute toxicity	
LD50 Oral - mouse - 2,314 mg/kg	
Skin corrosion/irritation	not available
Serious eye damage/eye irritation	not available
Respiratory or skin sensitization	not available
Chronic exposure	

No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC, ACGIH, NTP, or OSHA.

#### Potential health effects

Inhalation:	May be harmful if inhaled. May cause respiratory tract irritation.
Skin:	May be harmful if absorbed through skin. May cause skin irritation.
Eyes:	May cause eye irritation.
Ingestion:	May be harmful if swallowed.

### **Section 12 - Ecological Information**

Toxicity no data available Persistence and degradability no data available Bioaccumulative potential no data available Mobility in soil no data available PBT and vPvB assessment no data available Other adverse effects no data available

#### Section 13 - Disposal Considerations

ProductOffer surplus and non-recyclable solutions to a licensed disposal company.Contaminated packagingDispose of as unused product.

#### **Section 14- Transport Information**

DOT (US)	Not dangerous goods
IMDG	Not dangerous goods
IATA	Not dangerous goods

# **Section 15 - Regulatory Information**

OSHA Hazards Target Organ Effect, Irritant, Teratogen SARA 302 Components SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section

# **MSDS**

# **Material Safety Data Sheet**

Date prepared: May 28, 2012

### Date revised: N/A (first edition)

302.

#### SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

# SARA 311/312 Hazards

Acute Health Hazard, Chronic Health Hazard

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

#### Pennsylvania Right To Know Components

2-Naphthacenecarboxamide, 7-chloro-4-(dimethylamino)-1,4,4a,5,5a,6,11,12a-octahydro-3,6,10,12,12a -pentahydroxy-6-methyl-1,11-dioxo- CAS-No.64-72-2

#### New Jersey Right To Know Components

2-Naphthacenecarboxamide, 7-chloro-4-(dimethylamino)-1,4,4a,5,5a,6,11,12a-octahydro-3,6,10,12,12a -pentahydroxy-6-methyl-1,11-dioxo- CAS-No.64-72-2

#### California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

#### **Section 16 - Other Information**

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.