Material Safety Data Sheet

Date prepared: Sep 5, 2014 Date revised: N/A (first edition)

Section 1 - General Information

Chemical Name: Allethrin **Synonyms:** NA **CAS Reg. Number:** 584-79-2 **Catalog Number:** 1ST20484 Manufacturer's Name: A Chemtek, Inc. 100 Barber Avenue

Worcester, MA 01606, USA

Telephone Number: 508-471-4121 Fax Number: 508-845-9201

For R & D use only

Section 2 - Hazardous Ingredients/Identity Information

Classification of the substance or mixture

Acute toxicity, Oral (Category 4), H302 Acute toxicity, Inhalation (Category 4), H332 Acute aquatic toxicity (Category 1), H400 Chronic aquatic toxicity (Category 1), H410

Label elements



Pictogram Signal word

H302 + H332Harmful if swallowed or if inhaled

H410 Very toxic to aquatic life with long lasting effects.

Precautionary statement(s)

Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. P261

P264 Wash skin thoroughly after handling.

P270 Do not eat, drink or smoke when using this product. Use only outdoors or in a well-ventilated area. P271

P273 Avoid release to the environment.

P301 + P312IF SWALLOWED: Call a POISON CENTER or doctor/ physician if you feel

P304 + P340IF INHALED: Remove victim to fresh air and keep at rest in a position

comfortable for breathing.

P312 Call a POISON CENTER or doctor/physician if you feel unwell.

P330 Rinse mouth. P391 Collect spillage.

P501 Dispose of contents/ container to an approved waste disposal plant.

Section 3 - Composition/Information on Ingredients

Formula: C19H26O3 Molecular Weight: 302.41 g/mol 584-79-2 CAS-No.

Section 4 - First Aid Measures

If inhaled

If breathed in, move person into fresh air. If not breathing give artificial respiration.

In case of skin contact

Wash off with soap and plenty of water.

In case of eye contact

Flush eyes with water as a precaution.

Material Safety Data Sheet

Date prepared: Sep 5, 2014 Date revised: N/A (first edition)

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

Section 5 - Fire-Fighting Measures

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for fire-fighters

Wear self-contained breathing apparatus for firefighting if necessary.

Section 6 - Accidental Release Measures

Personal precautions Avoid dust formation

Environmental precautions Do not let product enter drains.

Methods for cleaning up Sweep up and shovel. Keep in suitable, closed containers for disposal.

Section 7 - Handling and Storage

Handling Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for

preventive fire protection.

Storage Keep container tightly closed in a dry and well-ventilated place.

Recommended storage temperature: 2-8 °C.

Section 8 - Exposure Controls/Personal Protection

Contains no substances with occupational exposure limit values

Personal protective equipment

Respiratory protection

For nuisance exposures use type P95 (US) or type P1 (EU EN 143) particle respirator. For higher level protection use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection For prolonged or repeated contact use protective gloves

Eye protection Safety glasses

Hygiene measures General industrial hygiene practice

Section 9 - Physical and Chemical Properties

Appearance:

Form liquid

Color not available

Safety data

pH: not available
Melting point: not available
Boiling point: not available
Flash point: not available
Lower explosion limit: not available
Upper explosion limit: not available
Water solubility: not available

Section 10 - Stability and Reactivity

Storage stability Stable under recommended storage conditions

Materials to avoid Strong oxidizing agents

Hazardous decomposition productsHazardous decomposition products formed under fire conditions.

- Carbon oxides.

Section 11 - Toxicological Information

Material Safety Data Sheet

Date prepared: Sep 5, 2014 Date revised: N/A (first edition)

Acute toxicity: LD50 Oral - rat - 685 mg/kg

LC50 Inhalation - mouse -> 2,000 mg/m3 LD50 Dermal - rabbit - 11,332 mg/kg

Irritation and corrosion:not availableSensitization:not availableRespiratory or skin sensitization:not available

Chronic exposure

IARC: 2B - Group 2B: Possibly carcinogenic to humans (Allethrin)

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a

known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a

carcinogen or potential carcinogen by OSHA.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a

carcinogen or potential carcinogen by OSHA.

Potential health effects

Inhalation: May be harmful if inhaled. May cause respiratory tract irritation. Skin: May be harmful if absorbed through skin. May cause skin irritation.

Eyes: May cause eye irritation.
Ingestion: May be harmful if swallowed.

Section 12 - Ecological Information

Elimination information (persistence and degradability)

Toxicity to fish LC50 - Cyprinus carpio (Carp) - 0.032 mg/l - 96.0 h Toxicity to daphnia and EC50 - Daphnia pulex (Water flea) - 0.021 mg/l - 48 h

other aquatic invertebrates

Ecotoxicity effects: not available **Further information on ecology:** not available

Section 13 - Disposal Considerations

Product Observe all federal, state, and local environmental regulations.

Contaminated packaging Dispose of as unused product.

Section 14- Transport Information

DOT (US)

Not dangerous goods

IMDG

UN number: 3082 Class: 9 Packing group: III EMS-No: F-A, S-F

Proper shipping name: ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S. (Allethrin)

Marine pollutant: Marine pollutant

IATA

UN number: 3082 Class: 9 Packing group: III

Proper shipping name: Environmentally hazardous substance, liquid, n.o.s. (Allethrin)

Section 15 - Regulatory Information

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

Acute Health Hazard, Chronic Health Hazard

Material Safety Data Sheet

Date prepared: Sep 5, 2014 Date revised: N/A (first edition)

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

Allethrin CAS-No. 584-79-2 Revision Date: 1989-12-01

New Jersey Right To Know Components

Allethrin CAS-No. 584-79-2 Revision Date: 1989-12-01

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16 - Other Information

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.

.