Material Safety Data Sheet

Date prepared: Sep 5, 2015 Date revised: N/A (first edition)

Section 1 - General Information

Chemical Name:FolpetCAS Reg. Number:133-07-3Catalog Number:1ST21234Manufacturer's Name:A Chemtek, Inc.

100 Barber Avenue

Worcester, MA 01606, USA

Telephone Number: 508-471-4121 **Fax Number:** 508-845-9201

For R & D use only

Section 2 - Hazardous Ingredients/Identity Information

Classification of the substance or mixture

Acute toxicity, Inhalation (Category 4), H332 Eye irritation (Category 2A), H319 Skin sensitisation (Category 1), H317 Carcinogenicity (Category 2), H351 Acute aquatic toxicity (Category 1), H400



Pictogram
Signal word
Warning

Hazard statement(s)

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

H332 Harmful if inhaled.

H351 Suspected of causing cancer. H400 Very toxic to aquatic life.

Precautionary statement(s)

P308 + P313

P201 Obtain special instructions before use.

P202 Do not handle until all safety precautions have been read and understood.

P261 Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.

P264 Wash skin thoroughly after handling.

P271 Use only outdoors or in a well-ventilated area.

P272 Contaminated work clothing should not be allowed out of the workplace.

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P304 + P340 + P312 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a

POISON CENTER or doctor/physician if you feel unwell.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove

contact lenses, if present and easy to do. Continue rinsing. IF exposed or concerned: Get medical advice/ attention. If skin irritation or rash occurs: Get medical advice/ attention.

P333 + P313 If skin irritation or rash occurs: Get medical advice/ attention.

P337 + P313 If eye irritation persists: Get medical advice/ attention.

P363 Wash contaminated clothing before reuse.

P391 Collect spillage. P405 Store locked up

P501 Dispose of contents/ container to an approved waste disposal plant.

Section 3 - Composition/Information on Ingredients

Formula: C9H4Cl3NO2S

Material Safety Data Sheet

Date prepared: Sep 5, 2015 Date revised: N/A (first edition)

Molecular Weight: 296.56g/mol CAS-No. 133-07-3

Section 4 - First Aid Measures

If inhaled

If breathed in, move person into fresh air. If not breathing give artificial respiration.

In case of skin contact

Wash off with soap and plenty of water.

In case of eye contact

Flush eyes with water as a precaution.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

Section 5 - Fire-Fighting Measures

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for fire-fighters

Wear self-contained breathing apparatus for firefighting if necessary.

Section 6 - Accidental Release Measures

Personal precautions Wear respiratory protection. Avoid breathing vapours, mist or gas. Ensure

adequate ventilation. Evacuate personnel to safe areas. For personal protection see section 8.

Environmental precautions Prevent further leakage or spillage if safe to do so. Do not let product enter

drains. Discharge into the environment must be avoided.

Methods for cleaning up Soak up with inert absorbent material and dispose of as hazardous waste. Keep

in suitable, closed containers for disposal.

Section 7 - Handling and Storage

Handling Avoid contact with skin and eyes. Avoid inhalation of vapour or mist.

Storage Keep container tightly closed in a dry and well-ventilated place.

Containers which are opened must be carefully resealed and kept upright to prevent leakage.

Recommended storage temperature: 2-8°C.

Section 8 - Exposure Controls/Personal Protection

Contains no substances with occupational exposure limit values

Personal protective equipment

Respiratory protection

Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multipurpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection For prolonged or repeated contact use protective gloves

Eye protection Safety glasses

Hygiene measures General industrial hygiene practice

Section 9 - Physical and Chemical Properties

Appearance:

Form solid

Color not available

Safety data

pH: not available Melting point: not available

Material Safety Data Sheet

Date prepared: Sep 5, 2015 Date revised: N/A (first edition)

Boiling point: not available
Flash point: not available
Lower explosion limit: not available
Upper explosion limit: not available
Water solubility: not available

Section 10 - Stability and Reactivity

Storage stabilityStable under recommended storage conditionsMaterials to avoidStrong oxidizing agents, Strong bases

Hazardous decomposition productsHazardous decomposition products formed under fire conditions.

- Carbon oxides, Nitrogen oxides (NOx), Sulphur oxides, Hydrogen chloride gas.

Section 11 - Toxicological Information

Acute toxicity:not availableIrritation and corrosion:not availableSensitization:not availableRespiratory or skin sensitization:not available

Chronic exposure

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

Potential health effects

Inhalation: no data available
Skin: no data available
Eyes: no data available
Ingestion: no data available

Section 12 - Ecological Information

Elimination information (persistence and degradability) not available

Ecotoxicity effects: not available **Further information on ecology:** not available

Section 13 - Disposal Considerations

Product Offer surplus and non-recyclable solutions to a licensed disposal company.

Contaminated packaging Dispose of as unused product.

Section 14- Transport Information

DOT (US):

Not dangerous goods

IMDG:

UN number: 3077 Class: 9 Packing group: III EMS-No: F-A, S-F

Proper shipping name: ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S. (Folpet (ISO))

Marine pollutant:yes

IATA:

UN number: 3077 Class: 9 Packing group: III

Material Safety Data Sheet

Date prepared: Sep 5, 2015 Date revised: N/A (first edition)

Proper shipping name: Environmentally hazardous substance, solid, n.o.s. (Folpet (ISO))

Section 15 - Regulatory Information

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

Folpet (ISO) CAS-No. 133-07-3 Revision Date

New Jersey Right To Know Components

Folpet (ISO) CAS-No. 133-07-3 Revision Date

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16 - Other Information

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.

.